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TESTING DECLARATION

Asal declare to test the dubnoff waterboath 750 with horizontal oscillation, c/o its seat in Cernusco S/N (Mi), before the sale.

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This manual is an integral part of the product with which it is enclosed and must be conserved for future consultation. Carefully read the recommendations and instructions contained in this instruction booklet before carrying out installation, operating, adjusting or maintenance operations, since they give important information regarding the SAFETY OF USE and MAINTENANCE.

Even if this device is provided with protections, the use of equipment without following the described procedures could produce a wrong working and equipment's damages.

CERTIFICATO

Nr. 50 100 3290 Rev.006

SI ATTESTA CHE / THIS IS TO CERTIFY THAT

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È CONFORME AI REQUISITI DELLA NORMA
HAS BEEN FOUND TO COMPLY WITH THE REQUIREMENTS OF

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Design, production and service of scientific appliance and laboratory equipment. Trade of laboratory equipments, disposable products, glassware and medical items for scientific and technical purposes (IAF 19, 29)

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Dal / From: 2021-08-02

Al / To: 2024-07-25

Andrea Coscia
Direttore Divisione Business Assurance
Business Assurance Division Manager

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Refer to Fig. 1

- 1 Main switch bright
- 2 Resistance ignition led
- 3 In tank liquid temperature display
- 4 Failure probe led
- 5 Low water level led
- 6 Temperature calibration button
- 7 Setted temperature display
- 8 Up-down temperature value regulation buttons
- 9 Programmable channels buttons
- 10 Cycle starting and stopping button
- 11 Stirrer speed controller
- 12 Stirrer timer

Fig. 2

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- 13 Discharge cock water
- 14 Discharge water
- 15 Heating overload fuses
- 16 Shaking overload fuses
- 17 Power cable

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
file: user guide
mod. 750
page: 18
release: 09/21

10 – DEACTIVATION AND DISMANTLING - Compliance with WEEE Directive

Italian legislative decree 14 marzo 2014, n. 49, implementation of WEEE Directive 2012/19/UE, rules the criteria for disposal and recovery of electrical and electronic materials. All the products sold by ASAL s.r.l., are professional equipments, belong to the category 8 annex 1 Dlgs 14 marzo 2014 n. 49, so they shall be subject to the directive. The product should not be disposed of as household waste not recyclable. Use the appropriate structures of local WEEE collection for disposal of this product or call dealer or manufacturer, subscribed to the “National Register of subjects required to finance the management systems WEEE” and the collective system for the management of WEEE called “Consorzio ECOLIGHT”.

The manufacturer is obliged to retire if the equipment produced after 13 August 2005; if produced before that date, the withdrawal is subject to the purchase by the customer for an equivalent product, and used the same functions (art. 24 Dlgs 14 marzo 2014 n. 49).

Disposal inappropriate may have negative consequences for the environment and human health. The recycling of materials helps conserve natural resources.

The disposal of WEEE as municipal waste is punished with sanctions specifically provided for in the case of illegal disposal of this waste.

This symbol on the product or packaging indicates that the same once finished, its use must be delivered to a collection centre for the recycling of electrical and electronic equipment.

Only trained and equipped personnel have to do these operations, and they must be in accordance with Dlgs 14 marzo 2014, n. 49 and European directive 2012/19/UE for the treatment and recycle of the electric and electronic wastes.

Disconnect the power rated (see the paragraph 3 - installation).

Strip these parts:

1. Electric parts
2. Plastic parts
3. Plexiglas parts
4. Iron parts
5. Syntetic parts

Follow the instructions of one's municipality.

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7 - GUARANTEE

- ◆ The good working of the described product is guaranteed from ASAL for a period of 12 months from the purchase date.
- ◆ The guarantee includes the repair or the free exchange of the instrument components that ASAL technicians admit to be defective in terms of manufacture or material (carriage and packing costs to be paid by the buyer).
- ◆ The guarantee will be invalidated by negligence, wrong installation and uses which are not observed in this instructions book, incident and tampering or transfer damage (excluding when it is made by our authorized people).
- ◆ The guarantee extension is no longer valid after repair.
- ◆ In case of restitution, the equipment has to be returned clean in order to avoid any problem for our personnel. We will refuse to draw back all the instruments that do not observe these sanitary regulations.
- ◆ Fuses are not included in the guarantee.

8 - DOCUMENT VALIDITY

ASAL reserves the right to review this document without any notice (any variations do not prejudice the observation of the directives).

9 - SPARE PARTS

For any spare parts needed contact the builder (see address and telephone number on the last page).

file: user guide
mod. 750
page: 16
edition: 09/21

5 - MAINTENANCE

To Clean the box powdered painted use domestic detergents or cleaners.
For the parts in stainless steel it is suggested use soapy water, detergents for glass or denatured alcohol, however non-abrasive products.

Do not use products with chlorine for parts in stainless steel.

The maintenance has to be carried out by experts.
To look at the good state of the item and plug.

In order to avoid any contact with hot parts, the operations have to be done in total safety conditions (disconnect the plug from the supply power cable).

6 - WORKING DEFECTS

In case of lack of tension, item stop working.

To work again follow the instructions:

- ◆ if the timer was positioned in "MANUAL CONTINUOUS", item automatically start working.
- ◆ if timer was switched on check that the time siglaed by the timer is in accordance with the operator needed.

Push START to heat, the old memorized temperature's values are in use.

In case of failure to ignite:

- ◆ to cut out the plug and check that the fuses are in a good state
- ◆ check the outlet is in a good state

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4.1 - INSTRUMENT CALIBRATION

The panel has a digital calibration with a max. adjustable range of +/- 2.0 °C, by 0.1 steps.

In order to make the instruments calibration see the below instructions:

- ◆ Set a channel on a temperature on which the maximum precision is required (Pos. 8, Fig. 1).
- ◆ Ignite the heating pressing the "START / STOP" key (Pos. 9, Fig. 1).
- ◆ Wait for the temperature stabilization.
- ◆ After 10 minutes, compare the temperature shown on the panel display with the thermometer before placed inside the tank: for example, if the calibration has been done at 40.0 °C (so the temperature displayed is 40.0 °C) and the thermometer shows 39.6 °C, it is necessary to set the ADJ value until the main display will show 39.6° C (Pos. 2, Fig. 1).
- ◆ Press the "ADJ" button (Pos. 6, Fig.1): the SET displays blink showing the back calibration value.
- ◆ Press slowly the up/down arrow key (Pos. 10, Fig.1) until the reference thermometer IDENTICAL value is reached.
- ◆ Press again the "ADJ" button (Pos. 6, Fig. 1).

4.2 - ERROR CODES

OP	PROBE IS OPEN	PT 100 circuit is open
CC	PROBE IN SHORT CIRCUIT	PT 100 circuit is in short circuit

file: user guide
mod. 750
page: 14
edition: 09/21

4 - INSTRUCTION FOR USE

WARNING: do not use distilled water, only deionized or descaled water, to prevent the limestone formation and incrustations.

The water level to be inserted depends on the size of containers used. Enter a quantity of water suitable for heating containers. Maximum level at about 5 cm. from the top edge of the tank.

Refer to Fig. 1 on page 20

Turn the main switch on, green led light on (Pos 1, Fig.1).

- STIRRER CONTROLS

If it is necessary to use the stirrer in a continuous way, put the timer on "Manual Continuous" (Pos.3, Fig.1).

If it is necessary to use the stirrer for a defined period of time, between 1 and 60 minutes, put the timer knob on "Timed" for the required period (Pos.3, Fig.1).

In order to settle the oscillation's speed, use the "Speed Controller" (Pos.2, Fig 1).

- WATERBATH CONTROLS

At the starting point the panel makes a general Lamp Test of 2 seconds and then it is ready to work. Choose the suitable working channel using the CH button (Pos. 8, Fig. 1).

The heating starts to work at the pre-set value pressing "STAR / STOP" (Pos. 9, Fig. 1). During the heating a LED shows the resistance ignition (Pos. 3, Fig. 1).

To stop the cycle press the "START / STOP" button (Pos. 9, Fig. 1).

To programme the visualized working channel, press the up/down arrow key (Pos.10, Fig.1) until the desired temperature will be reached. At the end of the programming phase the read in is automatically done.

To choose a channel press the "CH 1,2 or 3" button (Pos.8, Fig.1).

While the heating is working it is possible to make every kind of changes, both of the SET value and the channel. Each change is stored.

file: user guide
mod. 750
page: 11
edition: 09/21

TABLE OF CONTENTS

1	DESCRIPTION	page 12
2	OPERATING FEATURES	page 12
3	INSTALLATION	page 13
4	INSTRUCTION FOR USE	page 14
5	MAINTENANCE	page 16
6	WORKING DEFECTS	page 16
7	GUARANTEE	page 17
8	DOCUMENT VALIDITY	page 17
9	SPARE PARTS	page 17
10	DEACTIVATION AND DISMANTLING – Compliance with WEEE Directive 2012/19/EC D.lgs 14 marzo 2014, n. 49	page 18
11	WIRING DIAGRAM	page 19
12	PARTS LIST	page 20
13	ACCESSORIES	page 22
14	EU DECLARATION OF CONFORMITY	page 23
15	TABLE OF SYMBOLS AND TESTING DECLARATION	page 24

file: user guide
mod. 750
page: 12
edition: 09/21

1 - DESCRIPTION

Dubnoff waterbath with horizontal oscillation is widely used in the chemical, industrial and pharmaceutical fields. Dubnoff unit maintains test tubes, beakers and flasks in constant stirring

Condensate-free hinged cover lid in stainless steel Aisi 304. Cock for water drain. Outer casing made of powdered painted steel acids resistant. Internal tank in stainless steel Aisi 304.

Inner basket and work plane with horizontal oscillation, in stainless steel Aisi 304.

Main switch bright.

WARNING:

DO NOT USE DISTILLED WATER, ONLY DEIONIZED OR DESCALED WATER

2 – OPERATING FEATURES

2.1 - TECHNICAL SPECIFICATIONS

- ◆ Electronic regulation of oscillation in range 40 to 150 per minute. Oscillation amplitude: 14 or 24 mm.
- ◆ It is possible to utilise it in a continuous way or for a preset period up to 60 minutes by using the timer control.
- ◆ HCMOS technology microprocessor controlled.
- ◆ Liquid temperature in the tank measurement by PT100 thermal sensitive probe, visualisation on left display with resolution 0.1°C, precision +/- 1 °C .
- ◆ Set-point temperature indication on right display, settable from 10.0 to 99.9 °C by steps of 0.1 °C.
- ◆ Three programmable channels directly recalled by operator.
- ◆ Digital calibration of +/- 2 °C.
- ◆ Alarm indication of failure probe and low water level.
- ◆ Indication of heating resistance operation.
- ◆ Antisplash protection degree IP 54.
- ◆ Test panel following the EMI IEC 801.4 regulations.
- ◆ Discharged tap
- ◆ Stainless steel inner tank suitable to contain 2 test leads with 36 seats diam.mm. 18 or other containers
- ◆ Work area dimensions WxDxH: 360x260x140 mm. = Capacity 13 litres. Maximum height useful with closed cover 220 mm.
- ◆ Work surface dimensions WxD: 360x260 mm.
- ◆ Overall dimensions WxDxH: 770x400x380 mm.
- ◆ Weight: 20 Kg.

file: user guide
mod. 750
page: 13
edition: 09/21

2.2 - ELECTRICAL SPECIFICATIONS

- ◆ Power supply: 230 V - 50 Hz.
- ◆ Electrical input: 1530 W
- ◆ Overload fuses: 2 x 10 AF (5x20) mm. + 2 x 0.5 AF (5x20) mm.
- ◆ Connection outlet: 16 A

2.3 - CLIMATIC CONDITION FOR STORAGE, PACKING AND USE

- ◆ Altitude up to 2000 m
- ◆ Temperature 5 °C to 40 °C
- ◆ Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C
- ◆ Mains supply voltage fluctuations not to exceed ± 10% of nominal voltage
- ◆ Transient overvoltages II
- ◆ Pollution degree 3

3 - INSTALLATION

In order to have a functional working position, it is necessary to put the item on a stable support. To move it lift from lower part.

Allow at least 5 cm. between the water discharge located on the rear wall of the bathroom and the nearest wall.

Don't use with human health dangerous materials. Don't use in dangerous atmosphere. Don't use with flammable and dangerous materials or toxic gases that can bring about implosion or explosion, without precaution.

Asal s.r.l. is not aware of the substances contained in vessels or tubes in the heated water bath model 750.

You alone are responsible for the handling of these substances

Before the connection of the equipment check that:

- ◆ the rated voltage value indicated on the equipment nameplate corresponds to the actual mains values.
- ◆ the power supply has an effective connection on ground.
- ◆ the electric controls are in a no-working position.
- ◆ plug in the power supply cable.