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Noi

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Declaration of compliance

We

manufacturer VELP SCIENTIFICA s.r.l.
address Via Stazione, 16
USMATE (MI)
Italy

Under our responsibility declare that the product conforms to enforce prescriptions on safety and hygiene of machines directive 89/392/CEE, low voltage directive 73/23/CEE, electromagnetic compatibility directive 89/336/CEE plus modifications and that the documents listed in the annex V are available at Velp's offices as per the machine directive.

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cod. 10.0162 (220-240 V / 50 Hz)

Instruction for use: magnetic stirrer type ARE

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IMPORTANT:

READ THE INFORMATION CONTAINED IN THE PRESENT MANUAL BEFORE USING THE UNIT.

THE MANUFACTURER DOES NOT ACKNOWLEDGE ANY RESPONSIBILITY FOR AN IMPROPER USE OF THE EQUIPMENT, NOT RESPONDING TO DIRECTIONS FOR USE

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The labels applied to the unit warn the user on the dangers to which he is exposed during the use or the maintenance. The labels must be left on the unit and substituted if they are no more readable.

Warn of danger

To prevent accidents which could be caused by the unit

Start the operation

Do not start to use the unit before reading directions for use and maintenance

SAFETY REGULATIONS

- 1) The heating plate, if programmed, can reach a temperature of up to 370°C, not only during the heating phase but also in the cooling down phase when the appliance is switched off.
- 2) The magnetic fields are active on magnetic or metallic pieces (magnetic supports, CDs, floppy discs, cardiac stimulators, magnetic cards). Keep these pieces far away from the magnetic stirrer plate and the magnetic stir bars.
- 3) The stirrer must not be used:
 - in atmospheres at risk
 - in bain-marie
 - to stir combustible liquids that have a low combustion temperature
- 4) The container for the material being used must be smaller in diameter than that of the working surface to prevent any risk of tipping over.
- 5) If the heating of the heating plate is used the container of the product being used must be compatible with the temperature used. Normally, in these cases, Pyrex glass containers are used.

CLEANING

Always unplug the unit before cleaning. Cleaning must be done when the heating plate has cooled down, using a cloth dampened with non inflammable and non aggressive detergents.

PERSONNEL PROTECTION EQUIPMENT

The personnel protection equipment must be compatible with the temperatures reached, the risks deriving from the material being used and glass containers in compliance with the before mentioned safety regulations.

MAINTENANCE

In compliance with the law regarding product warranty, repairs of our units must be done on our premises, barring other agreements with local distributors.

WARRANTY

The warranty is effective from the date shown on the delivery note and refers to the registration number of the individual article.

The manufacturer, in an attempt at constantly improving the quality of his products, reserves the right to vary the characteristics without prior notice.

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1. Introduction

The magnetic stirrer is the basic unit of any laboratory.

The **ARE** model is designed for general use, therefore, for all those applications that require programming of the stirring speed and the temperature of the heating plate.

It is used to efficiently mix a liquid using the driving magnet of a stir bar positioned inside the container. The size and shape of the magnetic stir bar determines the efficiency of the stirring itself in line with the stirring speed. It is advisable, therefore, to choose the most suitable magnetic stir bar in relation to the quantity of liquid being used as well as the type of container used.

The stir bar which satisfies most applications is code Velp 11.1060 (10x40 mm).

The stirring speed is adjusted using an impulse microprocessor which makes it possible to obtain optimum speed progression.

The temperature of the heating plate is controlled by an analogical thermostat which, using the sensor bulb, allows for a constant control of the temperature.

N.B. If the heating of the heating plate is used the container of the product being used must be compatible with the temperature used. Normally, in these cases, Pyrex glass containers are used.

2. Description of the unit

The use of new technologies has led to the creation of modern laboratory magnetic stirrers whose basic characteristics are safety, reproducibility, high performance levels and energy saving.

The structure of the unit, made out of pressure die-cast aluminium, is treated with epoxy resins that are resistant to the manifold chemical aggressions found in laboratories.

The meticulously studied pressure die-casting line, combines the fluidity of shape with the ergonomics of the unit:

- The socket for moving the unit is assisted by the wide central recess. When the unit is not connected to the electricity supply it can be transported using the handle on the back which is moulded directly into the pressure die-casting.

- The control panel is easily accessible, far removed from any dangerous heat sources for the operator and protected from any possible leakages of the liquid by a special conveyance channel.
- The particular geometry of the pressure die-casting means that no accidental overflowing of the liquid being used can enter into the internal parts of the unit.
- The electrical connections are positioned inside on the back and provides optimum electrical protection against the overflowing of liquid into the outer electrical parts.
- The engineering of the unit allows for an easy internal inspection for any possible technical assistance interventions.

The heating plate, made out of cast aluminium, is treated with a particular process which provides:

- Optimum heat distribution with high specific power thanks to the circular configuration
- Optimum even temperature in all points on the plate
- High resistance to thermal shock
- High resistance to chemical aggressions
- High resistance to surface abrasions

N.B.: The use of the heating plate with high temperatures may cause variations of the superficial colour but this doesn't change the characteristics of thermal, mechanical and chemical resistance.

The magnetic stirring is generated by the Velp Scientifica Patented Composed Magnet. It is driven by an asynchronous mono-phase brushless motor in alternate current which allows for a practically unlimited duration.

The stirring speed is adjusted using an impulse microprocessor which makes it possible to obtain optimum speed progression.

The speeds can be varied using the special knob from between 50 and 1.200 rpm.

The temperature control of the heating plate control is done electronically by means of an analogical bulb thermoregulator fitted directly into the heating plate.

This makes it possible to obtain a constant control of the heating plate temperature and all the set temperatures.

3. Checking the unit on receipt

The unit is supplied packaged in a special anti-impact container, and comes complete with:
- supply connection cable

4. Assembling and installation

Check that the information on the plate on the unit corresponds to that of the electricity supply. Connect the unit to the electricity supply using the appropriate cable supplied with the unit. Turn on the unit using the general switch. When on it lights up.

5. Start up

After having turned on the general switch position a beaker with a magnetic stir bar. To select the stirring speed turn the corresponding knob (**Stirrer**). The stirring speeds can be programmed between 50 and 1.200 rpm. To select the temperature of the heating plate turn the corresponding knob (**Heating**). The temperature can be programmed between 50 and 370°C.

6. Operations at the end of work

When the unit is not be used it is advisable to turn off the general switch. When off the light goes out

7. Maintenance

The unit does not require any ordinary and/or extraordinary maintenance if not the cleaning as described in this manual in the general introductory notes.

It may be necessary to replace the fuses. The unit is fitted with two fuses, found in the socket positioned on the back. To replace one or more of the fuses remove the connection socket and, using a screwdriver, lift up the small cover on the fuse box.

8. Disposal

To dispose of and get rid of the material and components of the unit follow the instructions and laws in force in the country where the unit is to be disposed of. Particular attention must be paid to the permanent ferrite magnets which make up the driving magnet system found at the top of the driving shaft.

9. Accessories

Magnetic stir bar diam. 6x35mm	11.1056
Magnetic stir bar diam. 10x40 mm	11.1060
Hemispheric bowl for 250 ml flasks	11.1071
Hemispheric bowl for 500 ml flasks	11.1072
Hemispheric bowl for 1000 ml flasks	11.1073


10. Spare parts

Description	10.0162 10.0162-60	10.1162 10.1162-20
Heating plate	10.0162.20	10.1162.20
Medium composed magnet	10.0162.06	
Electronic speed control card	EESH163	EESH164
Motor	MOMO076	MOMO086
5x20 fuses	ELFU004	ELFU005
Thermostat	TRTE001	
Ceramic spacer H 7.5 mm	MERN009	
Ceramic spacer H 4.5mm	MERN010	

11. Specifications

Dimensions BxHxP	mm	165x115x280
Weight	Kg	2,7
Overall power	W	780
Power of the heating plate	W	750
Diameter of the heating plate	mm	155
Programmable temperature range	°C	50 – 370
Type of temperature control		Analogical
Stirring power	W	30
Stirring capacity	Litres of H ₂ O	15
Programmable speed range	rpm	50 ÷ 1200
Type of motor control		electronic

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13. Suggestions

Any suggests and advice would be greatly appreciated to enable us to constantly improve this instruction manual.

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Pompe

POMPE PERISTALTICHE

** distribuiti solo in Italia - distributed only in Italy*

Thanks for choosing a VELP's product !

Since 1983 Velp offers to the professional in the field a range of sophisticated and reliable equipment which make available high levels of know-how and operative capacity at competitive price. Velp works according to ISO 9001 Certification of Quality. Instruments are built according to the International norms IEC 1010-1 and to the rules of CE mark.

We present you our product Lines:

Environment

C.O.D.
B.O.D. – B.O.D. Sensor
REFRIGERATED THERMOSTATS
INCUBATORS
JAR TEST
METALS IN SLUDGES AND SEDIMENTS
NEPHELOMETRIC TURBIDITY
NITROGEN IN DIFFERENT FORMS
PERFORMANCE OF ANAEROBIC DIGESTERS
PHENOLS
LEACHABILITY TESTS

Stirring

MAGNETIC STIRRERS
HEATING MAGNETIC STIRRERS
MULTIPLACE MAGNETIC STIRRERS
LIGHTED MAGNETIC STIRRERS
HEATING PLATES
VORTEX
STIRRERS

Food and feed

DIGESTERS
STEAM DISTILLING UNITS
SCRUBBER
CATALYST TABLETS
NITROGEN-FREE WEIGHING BOATS
EXTRACTORS FOR RAW FIBER
EXTRACTORS FOR DIETARY FIBER
EXTRACTORS FOR SELECTIVE EXTRACTION OF SOLUBLE PRODUCTS BY SOLVENTS (FOR FAT)

Pumps

PERISTALTIC PUMPS

e-mail: service@velp.it