

Gilson Test Report

IDENTIFICATION	
Model:	Pipetman L P20L
Useful volume range:	2 µL to 20 µL
Serial N°:	TK71515
Tips:	Diamond D200
TEST CONDITIONS	
Balance n°:	88
Sensitivity:	10 <sup>-6</sup> g
Correction factor Z (µL/mg):	1.0032
Correction for evaporation:	0.018 µL
Basis of adjustment:	Ex
Air temperature:	21.5°C
Barometric pressure:	101 kPa
Relative humidity:	60 %
CONTROL	
Date: 6 SEP 2022	Inspected by: LM

VOLUMETRIC DATA	
Gravimetric results are converted to volumetric data using the corresponding correction factor Z and taking into account the evaporation if applicable.	
Mean:	2.004 µL
Systematic error:	0.004 µL
es %:	0.20 %
Random error:	0.012 µL
CV%:	0.60 %
TEST RESULTS	
Selected volume:	2 µL — PASS
Nominal volume:	20 µL — PASS
WARRANTY PERIOD	
3 years	

Gilson Warranty

Gilson warrants this pipette against defects in material under normal use and service for the period specified in the warranty document, from the date of purchase. This warranty shall not apply to pipettes which are subject to abnormal use and/or improper or inadequate maintenance (contrary to the recommendations given in the user's guide), including, but not limited to pipettes which have been subjected to physical damage, improper handling, or spillage or exposure to any corrosive environment. This warranty shall also be void in the event pipettes are altered or modified by any party other than Gilson or its designates. Gilson's sole liability under this warranty shall be limited to, at Gilson's sole option, repair or replacement of any defective components of pipettes or refund of the purchase price paid for such pipettes.

The foregoing warranty is exclusive and Gilson hereby disclaims all other warranties, whether express or implied, including any warranties of merchantability and any warranties of fitness for a particular purpose, under no circumstances shall Gilson be liable for any consequential, punitive, indirect or incidental damages arising out of any breach of any express or implied warranty.



LT803115/H - © June 2016



Warranty registration card

To ensure registration of warranty, please complete and return this card, promptly to your local Gilson distributor, or online: [www.gilson.com/register](http://www.gilson.com/register). PLEASE PRINT CLEARLY.

Model:	P20L
Serial Number:	TK71515
Date:	6 SEP 2022

LAST NAME		FIRST NAME	
COMPANY/INSTITUTION			
DEPARTMENT			
MAILING ADDRESS			
CITY		STATE/COUNTRY	ZIP/POSTAL CODE
( )		( )	
PHONE	FAX	E-MAIL	
DATE OF PURCHASE / /		DO YOU WANT TO RECEIVE NEW PRODUCT INFORMATION BY	
		E-MAIL ? <input type="checkbox"/> YES <input type="checkbox"/> NO	
		BY FAX ? <input type="checkbox"/> YES <input type="checkbox"/> NO	
BUSINESS TYPE (Select most appropriate)		<input type="checkbox"/> PHARMACEUTICAL <input type="checkbox"/> UNIVERSITY/MED SCHOOL <input type="checkbox"/> GOVERNMENT	
<input type="checkbox"/> INDUSTRIAL <input type="checkbox"/> RESEARCH INSTITUTE <input type="checkbox"/> FORENSICS		<input type="checkbox"/> CLINICAL CHEMISTRY/DIAGNOSTICS <input type="checkbox"/> ENVIRONMENTAL	
<input type="checkbox"/> BIOTECH <input type="checkbox"/> FOOD/BEVERAGE <input type="checkbox"/> OTHER			

For further details on the warranty, please contact your local Gilson distributor, or online: [www.gilson.com/register](http://www.gilson.com/register)



For your convenience, you can detach the bar code sticker.